



# Measurement of Specular Reflectivity




In specular reflectivity measurements, neutrons impinge upon the sample surface at an angle  $\theta_i$  and are scattered at an angle  $\theta_f$ . The incident and exit angles are equal ( $\theta = \theta_i = \theta_f$ ) and incremented together. The wave vector  $Q_z$  is defined as  $4 \pi \sin \theta / \lambda$ , where  $\lambda$  is the neutron wavelength. Above the critical angle  $\theta_c$  for total internal reflection, the data show finite-size fringes whose separation are inversely related to the film layer thickness. After subtraction of the off-specular background, these data can be fit (or inverted) to obtain a real-space profile of the scattering length density as a function of depth.


Data



Results



# Measurement of Off-Specular Reflectivity



Measurement of the off-specular reflectivity provide information about the length scale of in-plane structural correlations. For transverse- $Q_x$  scans (i.e., rocking curves), the scattering angle  $2\theta$  is held constant while  $\theta_i$  and  $\theta_f$  are varied equally in opposite directions ( $\theta_i + \theta_f = \text{const}$ ). Typically a narrow specular peak, evident at  $Q_x=0$ , can be separated from the underlying diffuse scattering which is broad. The width of the diffuse peak is indirectly related to the inverse of the coherence length  $\xi$  of the in-plane roughness.

