



# SEE YOU ON THE FLIPPED SIDE

The construction and testing of RF Coils used to flip  $^3\text{He}$  polarization




Joelle Baer



# Presentation Overview

- What is  $^3\text{He}$  used for at the NIST Center for Neutron Research (NCNR) ?
- How do you flip the  $^3\text{He}$  polarization?
- What are the steps in creating and testing such a device?
- Findings from this summer's work

# What is a $^3\text{He}$ cell?



- Spin dependent neutron absorption

K.P. Coulter et al, NIM A 288, 463 (1990)

# What is a $^3\text{He}$ cell?

- Back-filled with  $^3\text{He}$  and a small amount of  $\text{N}_2$
- Combination of distilled Rb/K
- Different cell characteristics




Cell for Typical Beam Size  
(Slider)



Wide Angle Cell (Reliance)

# What is polarized $^3\text{He}$ used for at the NCNR?





[www.ncnr.nist.gov/equipment/he3nsf](http://www.ncnr.nist.gov/equipment/he3nsf)



# Polarizing the Cell

## Spin Exchange Optical Pumping (SEOP)




- Near infrared laser
- Electrons are polarized
- Spin exchange with  $^3\text{He}$
- Lengthy process

# Polarization Analysis using $^3\text{He}$

Cross Sections:  $\uparrow\uparrow$   $\downarrow\downarrow$   $\uparrow\downarrow$   $\downarrow\uparrow$

Unconventional Set up:




- Probe magnetic properties from a sample
- Measure four cross sections
- Example:  $\downarrow\uparrow$

# Flipping the $^3\text{He}$ Polarization

- $B_0$  Field compensation on end caps
- Neutron shielding on N end
- Sleeve creates perpendicular  $B_1$  field
- No material in neutron beam



Vulcan





Hamilton




# Flipping the $^3\text{He}$ Polarization

- Oscillating B-field,  $B_1$  perpendicular to the static  $B_0$
- $^3\text{He}$  precess around  $B_0$
- Adiabatic Fast Passage (AFP) NMR


$B_0$  along  
cylindrical axis  
of solenoid



$$B_{eff} = \sqrt{(B_0 - \omega/\gamma)^2 + B_1^2}$$



# Why minimizing AFP loss is important?




To flip more often than 5 min (3.3 mHz), it is desired to have a flipping efficiency close to 0.9999 (loss of  $10^{-4}$ )


Efficiency determined by the Amplitude and Rate used in AFP NMR

# Building the RF Coil

- Sheet of Teflon rolled into a cylinder
- Hole for cell tip
- Slits in sine distribution



Sleeve in  
solenoid  
with  
constant  $B_1$   
due to slit  
distribution



# Building the RF Coil



- Wound with copper wire
- Down the inside, along the rim's edge, up the inside

# Building the FID NMR Coil

FID NMR Coil



- Free Induction Decay (FID) NMR
- Coil placed around cell tip
- Information about cell polarization and lifetime




# Testing the RF Coil

- Components to test:
  - Match the RF coil to the B field of the solenoids
    - For both 55kHz and 90kHz
  - Sweep amplitude for sleeve and solenoid pairs
  - Sweep rate for sleeve and solenoid pairs

# Testing the RF Coil

- Information about polarization from Amplitude (Amp) and Lifetime (T2)
- Take FID NMR, flip cell, take another FID NMR



# Testing the RF Coil

Trial #	# of flips	Amp Bef	Amp Aft	Phase Bef	Phase Aft	Aft/Bef	Loss	% Loss		Aft/Bef	Loss	%Loss
Rate: 650 kHz/s										Averages		
1	2000	.07615	.07154	220	224.9	0.93942	3.125E-05	0.00312		0.93485	3.4E-05	0.00337
2	2000	.07154	.0665	224.9	229.4	0.93027	3.614E-05	0.00361				

$$\text{Let } d = \frac{\text{Amplitude After}}{\text{Amplitude Before}}$$



$$\text{Loss} = 1 - \sqrt[\text{\# of flips}]{d}$$

Trial 2:

$$d = \frac{.0665}{.07154} = 0.93027$$

$$\text{Loss} = 1 - \sqrt[2000]{.93027} = .0000361$$

# Testing the RF Coil



# Parameter Results

- Weasley, Potter, and Granger all same size sleeve
- Variation in optimal settings
- Difference possible due to construction



OPTIMAL SLEEVE & SOLENOID PARAMETERS

Sleeve	Solenoid	Amp	Rate
Weasley	Gemini	2.7V	930kHz/s
Weasley	Pollux	2.9V	910kHz/s
Weasley	Vulcan	2.5V	-
		4.3 V	925kHz/s
Potter	Gemini	3.25 V	880kHz/s
Potter	Pollux	3.4 V	930kHz/s
Potter	Vulcan	2.5 V	980kHz/s
Granger	Gemini	3.6V	880kHz/s
Granger	Vulcan	3.7V	870kHz/s
Hagrid	Pollux	4.25V	650kHz/s



# Acknowledgements



Many thanks to  
The SURF Program  
SURF Directors  
NCNR Staff  
Wangchun Chen  
Shannon Watson  
Tom Gentile  
and Gordon Jones



# QUESTIONS